tomas®-basic set

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Instructions for use

Mode d'emploi

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Instructions for use tomas®-basic set

Dear Customer

Thank you for having chosen a quality product from Dentaurum.

In order to use this product at its best for you and your patients, it is important to study and follow these directions for use carefully.

The written instructions cannot cover all eventualities during operation. For this reason our Hotline is available to answer any other questions and ideas that may arise.

Due to a constant development we recommend, even when you use the same products frequently, to study the enclosed latest issue of the instructions for use or refer to our website at www.dentaurum.de.

1. Manufacturer

Dentaurum J. P. Winkelstroeter KG · Turnstraße 31 · 75228 Ispringen · Germany

2. With reference to our quality standards

Dentaurum ensures a faultless quality of the products manufactured by us. These recommendations are based upon our own experiences. The user himself is responsible for the processing of the products. Responsibility for failures cannot be taken, as we have no influence on the processing on site.

3. Basic instrument preparation

All **tomas®**-instruments including the **tomas®**-basic tray must all follow the same preparation procedure prior to initial use and after every application (disinfection, cleaning, rinsing, drying, care and maintenance, control, steam sterilization). The disinfection and cleaning of each **tomas®**-instrument and of the **tomas®**-basic tray should be carried out separately. The dry instruments are then placed in the respective position into the tray and the tray is to be closed. It is recommended to seal the tray. Then the steam sterilization is carried out.

Detailed information can be found in the **tomas®** instructions for use brochure (REF 989-631-20), as well as on the internet: www.dentaurum.de under „downloads“.

4. Explanations of symbols

	Gamma sterilized
	Non-sterile
	Lot-number (CH.-B)
	Expiration date...
	Single use only
	Please observe instructions for use
REF	Order-No.

If you have any questions, please contact the Dentaurum Orthodontic Hotline: +49 72 31 / 803-550.

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